

CO2-T

Fractional Ablative Laser

- Atrophic Vaginitis
- Vaginal Laxity
- Stress Urinary Incontinence
- Postpartum Perineal Trauma
- Cosmetic Laser Surgery
- Fractional and Traditional skin Resurfacing
- Skin Surgery
- Plastic and Aesthetic Surgery



1.Functions of Every Probes

Fractional Laser Probe



- 1.Skin resurfacing,skin rejuvenaion,wrinkle removal(superficial peeling)
- 2.Scar Removal:Acne scars,burnt scars ans surgery scars
- 3.Collagen Stimulation
- 4.Ablation and coagulation for the treatment of skin furrows and other textural irregularities
- 5.Dermatological Surgery(stretch marks)

Normal Laser Probe



- Flat warts removal
- Resection surgery
- Nevus/mole removal

Vaginal(VRL) Probes



- Contracting vaginal
- Labium whitening
- Solve the problem of dyspareunia and colpoxerosis

4. Technical Parameters

Parameters	
CO2 laser source	Mental RF laser emitter
Wavelength	10600nm
Power of laser emitter	30w
Pulsed radio frequency	0.530W
Screen	10.4 color touch LCD screen
Scan pattern size	0.1x0.1mm – 20x20mm
Spot size	0.05mm
Spot distance	0.1 -2.6mm adjustable
Pulse duration	0.1 -10ms adjustable
Lifetime of laser emitter	8-12 years
Cooling system	Air
Aiming light wavelength	650nm red semiconductor laser
Program language	English, Spain, Russian...nine languages
Voltage	110v/220v, 50/60Hz
New weight	15kg
Packaging Size	77*53*50cm

5. Interface, easy to operate

